

Family health care TRF-1000

Non-contact Infrared Thermometer



Measures instantly



Auto-Power off



Memory function



Sound ON/OFF



Object temperature measurement



Color LCD Display



Body temperature measurement



Non-contact Infrared Thermometer TRF-1000

Features

- Accurate and measures instantly
- 1 minute after showing the result, device is turned off automatically
- Automatic memory function, save up to 9 measured temperatures
- No probe filter need and easy to clean
- Ergonomic design for user
- LCD display shows normal, slight fever, high fever by color changes
Green: Normality (35.1°C ~ 37.4°C)
Orange: Slight fever (37.5°C ~ 38.5°C)
Red: High fever (> 38.5°C)
- FSTN LCD provides clear visibility from any different angle

How to use

• Body temperature measurement.

Hold the thermometer on 1~3cm away from the temple or center of the forehead between the eyebrows above, then press the measuring button. Keep the temperature sensor steady until the confirmation beep sound on.

* When the buzzer switch is set to off, press the measurement button and hold it for approximately 2 seconds.



Confirmation beep sound is that the body temperature has been measured, and the measurement shows on Color LCD Display. If the body temperature is too low or normal, a 'beep' will sound for approximately 2 seconds with green screen. The device will sound beep 4 times when the reading temperature is slight fever or high fever.

• Object temperature measurement

Keep a distance of 1-3cm between the thermometer and the object, press the measuring button. hold the position until the confirmation beep sound on after pressing the measuring button



When the buzzer is set to on, press the measurement button and hold it for approximately 2 seconds. Confirmation beep sound is that the body temperature has been measured, and the measurement shows on Color LCD Display.

Green: When measuring an object, green screen is turned on (0~100°C).

Specifications

Item	Content	Dimension
Product name/Model	WI-A Healthcare Non-Contact Infrared Thermometer TRF-1000	
Measurement method	Infrared detection /Body temperature measurement	
Basic feature	Body temperature Measurement. Object temperature measurement Memory Function Sound ON/OFF	
Display	LCD display	
Measuring range	Body temperature : 34°C ~42.2°C Object temperature : 1°C ~100°C	
Laboratory accuracy	Body temperature : 36°C ~39°C ± 0.2°C Other: ± 0.3°C Object temperature ± 2°C	
Display resolution	0.1°C	
Operation condition	15°C ~ 40°C / 15 ~ 95% RH (non condensing)	
Storag condition	-20°C ~ 55°C / 15 ~ 95% RH (non condensing)	
atmosphere pressure	700 - 1060hPa	
Software version	TRF-1000 Ver1.0.0	
IP	IP22	
Memory	Automatic memory function, save up to 9 results	
Power saving mode	Automatic power off after 60 seconds	
Power	DC 3V, Battery(LR03) X 2, at least 1,000	
Dimension	161.4mm x 54.3mm x 55.6mm	
Weight	93g	

